
#### FITTING AND APPLICATION FOR 2-POST COLLAR

AJUSTE Y APLICACIÓN DEL COLLARIN DE 2 BARRAS  
AJUSTEMENT ET MISE EN PLACE DE COLLIER À 2 TIGES SUPPORT  
POSIZIONAMENTO E APLICAZIONE DEL COLLARE CON 2 SUPPORTI  
ENCAIXE E APLICAÇÃO PARA COLAR DE 2 APOIOS  
PASSEN EN AANBRENGEN VAN KRAAG MET 2 SPIJLEN  
AVPASSNING OCH APLICERING AV NACKKRAGE MED 2 STÖDPUNKTER  
2 ÇUBUKLU BOYUNLUK İÇİN YERLEŞTİRME VE UYGULAMA  
ΕΔΑΡΜΟΤΗ ΚΑΙ ΠΥΘΜΙΣΗ ΤΟΥ ΚΟΛΛΑΡΟΥ 2 ΣΗΜΕΙΩΝ


#### SIZING

TALLA  
GRÖßENANPASSUNG  
TALLE  
MISURA  
TAMANHOS  
MATEN  
STORLEKAR  
BEDEN BELÛREME  
ΕΠΙΛΟΓΗ ΜΕΓΕΘΟΥΣ


#### VEST FITTING

AJUSTE  
AJUSTEMENT  
ANPASSUNG  
POSIZIONAMENTO  
ENCAIXE  
BEVESTIGEN  
AVPASSNING  
TAKMA  
ΕΔΑΡΜΟΤΗ

#### COLLAR FITTING

AJUSTE  
AJUSTEMENT  
ANPASSUNG  
POSIZIONAMENTO  
ENCAIXE  
BEVESTIGEN  
AVPASSNING  
TAKMA  
ΕΔΑΡΜΟΤΗ


#### FITTING AND APPLICATION FOR 4-POST COLLAR

AJUSTE Y APLICACIÓN DEL COLLARIN DE 4 BARRAS  
AJUSTEMENT ET MISE EN PLACE DE COLLIER À 4 TIGES SUPPORT  
POSIZIONAMENTO E APLICAZIONE DEL COLLARE CON 4 SUPPORTI  
ENCAIXE E APLICAÇÃO PARA COLAR DE 4 APOIOS  
PASSEN EN AANBRENGEN VAN KRAAG MET 4 SPIJLEN  
AVPASSNING OCH APLICERING AV NACKKRAGE MED 4 STÖDPUNKTER  
4 ÇUBUKLU BOYUNLUK İÇİN YERLEŞTİRME VE UYGULAMA  
ΕΔΑΡΜΟΤΗ ΚΑΙ ΠΥΘΜΙΣΗ ΤΟΥ ΚΟΛΛΑΡΟΥ 4 ΣΗΜΕΙΩΝ

#### PREPARATION

PREPARACIÓN  
PRÉPARATION  
VORBEREITUNG  
PREPARAZIONE  
PREPARAÇÃO  
VOORBEREIDING  
FÖRBEREDELSE  
HAZIRLIK  
ΠΡΟΕΤΟΙΜΑΣΙΑ


#### COLLAR FITTING

AJUSTE  
AJUSTEMENT  
ANPASSUNG  
POSIZIONAMENTO  
ENCAIXE  
BEVESTIGEN  
AVPASSNING  
TAKMA  
ΕΔΑΡΜΟΤΗ



#### VEST FITTING

AJUSTE  
AJUSTEMENT  
ANPASSUNG  
POSIZIONAMENTO  
ENCAIXE  
BEVESTIGEN  
AVPASSNING  
TAKMA  
ΕΔΑΡΜΟΤΗ

#### ADJUSTMENT

AJUSTE  
AJUSTEMENT  
REGULIERUNG  
REGOLAZIONE  
AJUSTE  
AFSTELLEN  
JUSTERING  
AYAR  
ΠΡΟΣΑΡΜΟΤΗ

Manufactured by:  
DeRoyal Industries, Inc.  
200 DeBak Lane  
Powell, TN 37849 USA  
888.938.7228 or (901) 865.938.7228  
www.deroyal.com

EMERGO EUROPE  
Prinsesgracht 20  
5114 AP The Hague  
The Netherlands

Part #0-18086  
Revised 2/2021

©2021 DeRoyal Industries, Inc. All Rights Reserved. DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are registered trademarks of DeRoyal Industries, Inc.

### CTO - 2-POST AND 4-POST CERVICAL THORACIC ORTHOSIS

SINGLE PATIENT USE	
	<b>NON-STERILE</b>
	<b>NOT MADE WITH NATURAL RUBBER LATEX</b>
<b>RX ONLY</b>	<b>FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.</b>

#### IMPORTANT INFORMATION

Please read all warning and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury or re-injury inherent with the use of any brace.

#### INTENDED USE

The CTO Cervical Thoracic Orthosis is intended to immobilize the cervical and upper thoracic spine.

#### WARNING

- A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber. Use only as directed.
- Prolonged usage can result in muscle atrophy in the user. The user may require (1) a weaning period from the brace to build muscle strength (2) regular physical therapy exercises to maintain the strength of the leg (3) a limitation on duration of use, or (4) some combination of these.
- Inspect brace for damaged or missing components before use. Discontinue use and consult your prescriber if the brace or its components break or become damaged.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.** Check regularly to ensure circulation is not compromised.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace. Discontinue use if pressure injuries develop.
- While not an absolute contraindication to the use of this brace, take special care if the user is diabetic or has poor circulation. These users may have decreased skin sensitivity and are at greater risk for poor peripheral circulation and pressure injuries.

#### INSTRUCTIONS FOR USE

##### FITTING AND APPLICATION FOR 2-POST COLLAR PREPARATION

- The 2-post collar provides support and motion restriction without discomfort from the posterior support posts.
- To prepare the 2-post CTO for fitting, remove the posted collar

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back and replace it with the unposted collar back. Replace the adjustment clamp.

- Remove the posted collar front from the vest front and undo the shoulder straps.

#### SIZING

- The patient should be lying flat on a bed with no pillows.
- Arms should be placed along side of the patient's body.
- Use XTW sizing guide to select the correct collar size.
- To size the thoracic belt, measure the chest circumference approximately 1" below the tip of the xiphoid process. If the circumference is less than 39", cut the belt at the notch between 41 and 44.
- Attach the belt by aligning the appropriate keyhole over the large aluminum rivets.
- Put firmly to seat the belt for a secure fit on the rivets. Repeat this on the opposite side.

#### VEST FITTING

- While maintaining appropriate spinal precautions, place back vest components next to the patient and guide them under the patient while depressing the foam mattress cover. Repeat to position the collar back under the patient's head. **NOTE:** If log rolling is required, have appropriate licensed practitioner perform this maneuver.
- Position the vest front by centering nipple line in the cutouts of the panel.
- Unattach the loop strap on each side and attach the thoracic belt. If it is necessary to change the belt circumference, move the belt symmetrically to the next size keyhole slot and secure.
- When the fit is snug, but not uncomfortable, lock in the position by applying the overlapping loop strap.
- Snap the shoulder straps into place but do not tighten.

#### COLLAR FITTING AND ADJUSTMENT

- Patient must be in neutral alignment for appropriate collar fit. Position the collar front on the patient and center.
- Position the back portion of the collar center on the patient's neck. With the front portion properly secured, use the bilateral hook and loop straps and fasten the portions together.
- Carefully tighten the collar until the collar flexes to contour to the patient. Minor adjustments can be made to the straps and front and back flex tabs to achieve desired comfort. The patient's chin should be flush or nearly flush with the end of the chin piece.
- Slide the support posts on the front of the CTO upward until the patient's head is properly positioned and supported. Tighten the screws on the adjustment clamp with the tool provided.
- Ensure that the collar and the vest panels fit symmetrically. Gently tighten the shoulder straps to remove any slack. Adjustments can be made using the adjustment clamp and shoulder straps to obtain an optimal fit for the patient. **NOTE:** Do not overtighten shoulder straps causing patient discomfort.

#### FITTING AND APPLICATION FOR 4-POST COLLAR PREPARATION

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- Use XTW sizing guide to select the correct collar size.
- To size the thoracic belt, measure the chest circumference approximately 1" below the tip of the xiphoid process. If the circumference is less than 39", cut the belt at the notch between 41 and 44.
- Attach the belt by aligning the appropriate keyhole over the large aluminum rivets.
- Put firmly to seat the belt for a secure fit on the rivets. Repeat this on the opposite side.

#### FOR UPRIGHT APPLICATION - COLLAR FITTING

- Loosely fit the collar on the patient. **NOTE:** Patient must be in neutral alignment.
- Make sure collar is properly centered and then tighten the collar with the bilateral hook and loop straps. Be sure to tighten it enough to ensure that the back panel flexes to fit the patient and all of the slack is removed. Minor adjustments can be made to the straps and front and back flex tabs to achieve desired comfort.

#### VEST FITTING

- Position vest front by centering the nipple line in the cutouts of the panel.
- Unattach the loop strap on each side and attach the belt. If it is necessary to change the belt circumference, move the rivet to the next size keyhole slot.
- When the fit is snug, but not uncomfortable, lock in the position by applying the overlapping loop strap.
- Snap the shoulder straps into place but do not tighten.

#### ADJUSTMENT

- Slide the support posts on the front of the CTO upward until the patient's head is properly positioned and supported.
- Tighten the screws on the adjustment clamp with the tool provided.
- Adjust the collar back in the same manner. Be aware that the higher the placement of the collar back, the more effective the motion restriction. Tighten the screws on the adjustment clamp.
- Ensure that both collar and vest panels fit symmetrically. Gently tighten the shoulder straps to remove any slack. Adjustments can be made using the adjustment clamp and shoulder straps to obtain an optimal fit for the patient. **NOTE:** Do not overtighten shoulder straps causing patient discomfort.

#### STORAGE AND TRANSPORT CONDITIONS

KEEP DRY

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

#### WARRANTY

DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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