EN: POLYDERM™ WOUND DRESSINGS

ES: APÓSITOS PARA HERIDAS POLYDERM™

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ENGLISH: POLYDERM™ WOUND DRESSINGS

	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

DESCRIPTION

The DeRoyal® Polyderm™ wound dressings are sterile, absorbent, polyurethane foam wound dressings. They are available with a breathable, clear, adhesive film border; a water resistant, acrylic border; or without any adhesive or border.

INTENDED USE

Polyderm™ wound dressings absorb wound drainage, promote a moist environment helpful to healing, and serve as a barrier against water and bacteria when used as a secondary dressing.

INDICATIONS

- Polyderm™ wound dressings are indicated for minor abrasions, minor lacerations, and minor scalds and burns.
- When under the care or direction of a licensed practitioner they are indicated for:
 - Dermat ulcers (e.g., pressure ulcers, leg ulcers, and other secreting lesions)
 - Second degree burns
 - Donor sites
 - Tube sites
 - More severe lacerations, cuts and abrasions
 - Other superficial wounds

CONTRAINdicATIONS

- Polyderm™ wound dressings are contraindicated for:
 - Third-degree burns
 - Heavily bleeding wounds
 - Individuals who are sensitive to or have had an allergic reaction to the dressing or its materials
 - Use during surgical procedures as a surgical sponge

WARNINGS

- Read and understand all warnings, cautions, and directions completely and carefully before use. Failure

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All Rights Reserved. DeRoyal, the DeRoyal logo, Improving Care. Improving Business., and Polyderm are trademarks or registered trademarks of DeRoyal Industries, Inc.