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Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

Part#74-14905
Revised 11/2020

©2020 DeRoyal Industries, Inc.

All Rights Reserved.

DeRoyal, the DeRoyal logo and Improving Care. Improving Business are trademarks or registered trademarks of DeRoyal Industries, Inc.

ESOPHAGEAL STETHOSCOPE WITH TEMPERATURE SENSOR

	DO NOT RESTERILIZE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE

The DeRoyal® Esophageal Stethoscope with Temperature Sensor is to be used for routine monitoring of core body temperature, as well as to monitor heart and respiratory sounds in an anesthetized patient by inserting the stethoscope tube into the esophagus.

PRODUCT DESCRIPTION

The stethoscope portion of the device consists of a small diameter plastic tube upon which a flexible blue cuff is placed. When the device is inserted into the esophagus, the cuff facilitates transmission of sound to a male luer fitting at the proximal end of the tubing to which standard stethoscopes may be attached. The stethoscope is available in 9, 12, 18, and 24 French sizes.

Secured within the lumen of the tubing is a temperature monitoring sensor. The connector at the proximal end interfaces with DeRoyal-branded cables for connection to a monitoring unit to display temperature. The temperature monitoring sensor is available in 400 or 700 series configurations for connection to a compatible monitoring unit. Consult the user manual for the patient monitor to ensure compatibility.

CONTRAINDICATIONS

The use of an esophageal stethoscope may be

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contraindicado en pacientes con esophageal diverticulum or stenosis, in neonates or small infants undergoing neck surgery, in patients undergoing a tracheostomy or insertion of an internal jugular catheter, and in patients with conditions deemed unsafe for insertion of the stethoscope into the esophagus.

! WARNINGS

- Tracheal or bronchial insertion of stethoscope may cause airway obstruction and may be fatal if uncorrected. Correct placement may be verified with a laryngoscope.
- DO NOT** use if package is open or damaged.

ADVERSE REACTIONS

While adverse reactions are rare, the following have been associated with use of an esophageal stethoscope with temperature sensor during insertion or while the device was still in use:

- airway obstruction
- aspiration pneumonitis
- bronchial insertion
- electrical burns and shock
- epitaxis
- esophageal abrasion
- esophageal perforation
- tracheal insertion
- trauma to pharynx

CAUTIONS

- If incorrect temperature readings or excessive response time is noticed, discard probe and replace.
- Intubate patient with tracheal tube before inserting the esophageal stethoscope.
- Inform the surgeon that an esophageal stethoscope has been inserted.
- Lubricate the stethoscope before insertion and use accepted medical techniques in the insertion and removal of the stethoscope.
- Take care to ensure that the cable and connector do not get wet and proper techniques are used during electrosurgical procedures to reduce radio frequency interference current and potential burn to the patient.
- An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.
- Do not intertwine the cables, especially the monitor cables with the electrosurgical unit's

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cables.

- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- If unusual or inconsistent temperature readings are noticed, check connections. If issue is unresolved, replace temperature probes.
- Check the temperature monitoring equipment user manual for information regarding compliance with Clause 5 of IEC 60601-1-2 – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests.
- The probe is intended to be used with a DeRoyal-branded interface cable. Use with an incompatible cable may affect performance.

DIRECTIONS FOR USE

- Remove the Esophageal Stethoscope with Temperature Sensor from its sterile package.
- Verify the compatibility of the probe, interface cable, and patient monitor before use.
- Lubricate the cuff with suitable water soluble lubricant.
- Insert the stethoscope and guide it into the esophagus.
- Listen for airway leaks when inserting. Invasion of the trachea may permit air to escape around the tracheal tube cuff. Do not seal the leak. Remove the esophageal stethoscope and reinsert. To verify proper placement, use a laryngoscope.
- Adjust the stethoscope's depth in the esophagus to maximize the heart and respiratory sounds.
- Connect the temperature probe to the appropriate reusable temperature cable from DeRoyal®. Connect the cable to the patient monitor. Secure the cable using the drape clip.
- Follow the directions for use of the patient temperature monitor instrument.
- Remove the probe prior to extubation of the tracheal tube following accepted medical technique.
- Disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull cable or wire.
- Dispose of probe according to facility protocol.

SAFETY IN MRI NOT EVALUATED

The Esophageal Stethoscope with Temperature

Sensor's safety in MRI has not been evaluated, hence, the probe is considered MR Unsafe.

OPERATING SPECIFICATIONS

Rated output range: 25 to 45° Celsius
The reference body site is the core body temperature.
The measuring site is the esophagus.
Accurate to +/- 0.2° Celsius

Time response

- The temperature probe has a direct mode of operation.
- The maximum heating transient time is 2 minutes 45 seconds.
- The maximum cooling transient time is 1 minute 45 seconds.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	DO NOT USE IF PACKAGE IS DAMAGED
	RELATIVE HUMIDITY: 85% NON-CONDENSING
	STORAGE TEMPERATURE: -25°C TO +55°C
	OPERATING TEMPERATURE: 25°C TO 45°C

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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