
DeRoyal uses the symbols and meanings from standards ISO 15523-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU), a complete symbols glossary is available online at deroyal.com/symbols or by contacting customer service.

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PT: SONDA DE TEMPERATURA NASOFARÍNGEA

Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

PART#74-14908
Revised 2/2025
©2025 DeRoyal Industries, Inc.
All Rights Reserved. DeRoyal and the DeRoyal logo are registered trademarks of DeRoyal Industries, Inc.

NASOPHARYNGEAL TEMPERATURE PROBE

	MR UNSAFE
	DO NOT RESTERILIZE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT
Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE
The DeRoyal® Nasopharyngeal Temperature Probe is to be used for routine monitoring of core body temperature.

INTENDED USERS
The DeRoyal® Nasopharyngeal Temperature Probe is intended to be used only by a licensed practitioner or trained healthcare professionals.

PATIENT TARGET GROUP(S)
The intended patient population is pediatric and adults.

PRODUCT DESCRIPTION
The DeRoyal Nasopharyngeal temperature probe is inserted nasally to monitor a patient's core body temperature. It is sterile, individually packed, and comes in 9 French size.

The temperature monitoring sensor is located at the tip of the probe and secured inside the lumen of the tubing. The Nasopharyngeal temperature probe is available with a soft sheath length of 3' and 18". It is available in YSI 400 series configuration and designed to be compatible with the patient monitor using the appropriate DeRoyal reusable interface cable. Consult the user manual for the patient monitor to ensure compatibility.

INDICATIONS
The DeRoyal® Nasopharyngeal Temperature Probe is indicated for the monitoring of temperature via the nasopharynx cavity.

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CONTRAINDICACIONES

- The use of this device in the nasopharynx cavity may be contraindicated in patients with esophageal diverticulum or stenosis, in neonates or small infants undergoing neck surgery and in patients undergoing a tracheotomy or insertion of an internal jugular catheter.
- Nasopharyngeal insertion is not recommended for patients who are taking anticoagulants or when manipulation of mucosa is undesirable.
- DO NOT** use rectally.

! WARNINGS

- To ensure accurate nasopharyngeal temperature readings, the Nasopharyngeal probe should be placed in contact with mucosa in the posterior nasopharyngeal cavity.
 - Inspect packaging for any damage and product for any damage, contaminate or missing parts prior to use.
- Use in the nasopharynx may cause epistaxis (nosebleed) which can have a higher incidence rate in pregnant patients, children with large adenoids, and patients with clotting disorders.
- This is a single use device. It is not intended for disinfection and/or subsequent re-use. Re-use could result in microbial contamination and potential patient injury.
- This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure.
- The mucosal temperature can vary based on exact location (i.e., anterior versus posterior mucosa) and respiratory cycle (i.e., inspiration versus expiration). The exact location and timing of the respiratory cycle should be noted and taken into consideration when using a nasopharyngeal temperature probe.

ADVERSE REACTIONS

While adverse reactions are rare, the following have been reported to be associated with the use of temperature probes during insertion or while the device was in use:

- airway obstruction
- aspiration pneumonitis
- bronchial insertion
- electrical burns
- epistaxis
- esophageal perforation
- tracheal insertion

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- trauma to pharynx
- mucosal abrasion
- infection/sepsis
- allergic reaction to materials

CAUTIONS

- Lubricate the probe with a medical grade water-soluble lubricant before insertion and use accepted medical techniques during insertion and removal of the probe.
- Improper use or positioning of the device can produce inaccurate temperature readings that could result in patient hyperthermia, hypothermia, or an inaccurate diagnosis.
- Take care to ensure that the cable and connector do not get wet and that proper techniques are used during electrosurgical procedures to reduce radio frequency interference current and potential burns to the patient.
- An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.
- DO NOT** intertwine the cables, especially the monitor cables, with the electrosurgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- If probe has inaccurate, unstable, or no temperature readings, discard and replace.
- The nasopharyngeal probe is intended to be used with DeRoyal-approved interface cables. Use with an incompatible cable may affect performance.

! SAFETY IN MRI NOT EVALUATED

The Nasopharyngeal Probe's safety in MRI has not been evaluated, hence, the probe is considered MR Unsafe.

DIRECTIONS FOR USE

- Remove the Nasopharyngeal Temperature Probe from its sterile package using aseptic technique.
- Verify the compatibility of the probe, interface cable, and patient monitor.
- Lubricate the probe tip with a suitable water soluble lubricant.
- Insert the probe and guide it through the nasal cavity into the nasopharynx cavity.
- Connect the temperature probe to the appropriate reusable interface cable from DeRoyal. Connect the cable to the patient monitor. Secure the cable using drape clip.

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- Follow the directions for use of the patient temperature monitor instrument.
- When temperature monitoring is completed, carefully remove the probe.
- Once temperature monitoring is no longer needed, disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire. **DO NOT** reuse.

OPERATING SPECIFICATIONS

- Operating Temperature: 25° C to 45° C
- Rated Output Range: 35° C to 42° C
- Rated Extended Output Range: 25° C to 35° C and 42° C to 45° C
- The reference body site is the core body temperature.
- The measuring site is the nasopharynx.
- Accuracy in Rated Output Range: +/- 0.2° C
- Accuracy in Rated Extended Output Range: +/- 0.3° C
- Time response:
 - The temperature probe has a direct mode of operation.
 - The maximum heating transient time is 40 seconds.
 - The maximum cooling transient time is 30 seconds.

DISPOSAL

Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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