DeRoyal uses the symbols and meanings from standards ISO 15523-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU). A complete symbols glossary is available online at deroyal.com/symbols or by contacting customer service.

NASOPHARYNGEAL TEMPERATURE PROBE

| | |
|--|--|
| | DO NOT RESTERILIZE |
| | DO NOT REUSE |
| | MEDICAL DEVICE |
| | STERILIZED WITH ETHYLENE OXIDE |
| | NOT MADE WITH NATURAL RUBBER LATEX |
| | FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER. |

IMPORTANT
Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE
The DeRoyal® Nasopharyngeal Temperature Probe is to be used for routine monitoring of core body temperature.

PRODUCT DESCRIPTION
Probe is inserted nasally. It is sterile, individually packaged, and comes in different sensor types. It has a low-friction surface finish to improve ease of insertion. The temperature sensor beneath the distal tip is highly accurate and designed to be compatible with the temperature monitor using the appropriate DeRoyal® reusable interface cable.

CONTRAINDICATIONS
• The use of this device in the nasopharynx cavity may be contraindicated in patients with esophageal diverticulum or stenosis, in neoplas or small infants undergoing neck surgery and in patients undergoing a tracheotomy or insertion of an internal jugular catheter.
• Nasopharyngeal insertion is not recommended for patients who are taking anticoagulants or when manipulation of mucosa is undesirable.
• DO NOT use rectally.

WARNINGS
• Not tested for use in MR environments. Use in an MR environment may cause burn to patient or other serious injury.
• Inspect packaging for any damage and product for any damage, contamination or missing parts prior to use.
• While adverse events are rare, the following have been reported to be associated with the use of temperature probes during insertion while the device was in use:
• airway obstruction
• aspiration pneumonia
• bronchial insertion
• electrical burns
• epistaxis
• esophageal abrasion
• esophageal perforation
• tracheal insertion
• trauma to pharynx

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- May cause epistaxis (nosebleed). In the pregnant patient, children with large adenoids, and patients with clotting disorders, the incidence of epistaxis is high.
- The mucosal temperature can vary based on exact location (i.e., anterior versus posterior mucosa) and respiratory cycle (i.e., inspiration versus expiration). The exact location and timing of the respiratory cycle should be noted and taken into consideration when using a nasopharyngeal temperature probe.

CAUTIONS
• Lubricate the probe with a medical grade water-soluble lubricant before insertion and use accepted medical techniques during insertion and removal of the probe.

• Take care to ensure that the cable and connector do not get wet and that proper techniques are used during electro-surgical procedures to reduce radio frequency interference current and potential risk to the patient.

• An adequate electro-surgical dispersive ground electrode close to the active surgical site should be properly connected.

• DO NOT intertwine the cables, especially the monitor cables, with the electro-surgical unit's cables.

• The operation of the patient temperature monitor may be temporarily affected during electro-surgical activations. Unusual temperature readings should be checked.

• If probe has inaccurate, unstable, or no temperature readings, discard and replace.

• The nasopharyngeal probe is intended to be used with DeRoyal-approved interface cables. Use with an incompatible cable may affect performance.

MR COMPATIBILITY

The DeRoyal® Temperature Probe has **NOT** been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Probe in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

DIRECTIONS FOR USE

1. Remove the Nasopharyngeal Temperature Probe from its sterile package.
2. Verify the compatibility of the probe, interface cable, and patient monitor.
3. Lubricate the probe tip with a suitable water soluble lubricant.
4. Insert the probe and guide it through the nasal cavity into the nasopharynx cavity.
5. Connect the temperature probe to the appropriate reusable temperature cable from DeRoyal. Connect the cable to the patient monitor. Secure the cable using drape clip.
6. Follow the directions for use of the patient temperature monitor instrument.
7. Disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire.
8. This probe is **NOT** reusable.
9. Discard according to facility protocol.

OPERATING SPECIFICATIONS

- Rated output range: 25° to 45° Celsius
- The reference body site is the core body temperature.
- The measuring site is the nasopharynx.
- Accurate to +/- 0.2° Celsius
- Time response:

• The temperature probe has a direct mode of operation.

STORAGE AND TRANSPORT CONDITIONS

| | |
|--|---------------------------------------|
| | KEEP DRY |
| | KEEP AWAY FROM SUNLIGHT |
| | DO NOT USE IF PACKAGE IS DAMAGED |
| | RELATIVE HUMIDITY: 85% NON-CONDENSING |
| | STORAGE TEMPERATURE: -25°C TO +55°C |
| | OPERATING TEMPERATURE: 25°C TO 45°C |

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

NOTES

• Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

• The DeRoyal® Nasopharyngeal Temperature Probe is to be used for routine monitoring of core body temperature.

• PRODUCT DESCRIPTION

Probe is inserted nasally. It is sterile, individually packaged, and comes in different sensor types. It has a low-friction surface finish to improve ease of insertion. The temperature sensor beneath the distal tip is highly accurate and designed to be compatible with the temperature monitor using the appropriate DeRoyal® reusable interface cable.

• CONTRAINDICATIONS

• The use of this device in the nasopharynx cavity may be contraindicated in patients with esophageal diverticulum or stenosis, in neoplas or small infants undergoing neck surgery and in patients undergoing a tracheotomy or insertion of an internal jugular catheter.

• Nasopharyngeal insertion is not recommended for patients who are taking anticoagulants or when manipulation of mucosa is undesirable.

• DO NOT use rectally.

WARNINGS

• Not tested for use in MR environments. Use in an MR environment may cause burn to patient or other serious injury.

• Inspect packaging for any damage and product for any damage, contamination or missing parts prior to use.

• While adverse events are rare, the following have been reported to be associated with the use of temperature probes during insertion while the device was in use:

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