

WARRIOR® WRIST LACER

	NOT MADE WITH NATURAL RUBBER LATEX
	SINGLE PATIENT USE ONLY
	NON-STERILE
RX ONLY	Federal U.S.A. law restricts this device to sale or use by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please read all instructions, warnings and precautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE

The DeRoyal® Warrior® Wrist Lacer brace is intended to provide compression, support and alignment to the wrist

INDICATIONS





Carpal tunnel syndrome; Tendonitis; Sprains and strains of the wrist; Overuse injuries

CAUTION

- This product is to be fitted initially by a physician or properly licensed practitioner who is familiar with the purpose for which it is intended. The physician or practitioner is responsible for providing wearing instructions and precautions to other healthcare practitioners or healthcare providers involved in the patient's care as well as the patient. Consult your physician or properly licensed practitioner immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this product. Discontinue use if pressure ulcers develop. Follow the instructions of your physician or properly licensed practitioner for duration of use. These instructions do not supersede hospital protocol and/or direct orders of patient's physician.
- Discontinue use of the device and consult physician or properly licensed practitioner if it becomes damaged in any way, limiting its ability to perform as intended.
- Inspect brace for damaged or missing components prior to patient application.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.**

INSTRUCTIONS FOR USE

1. Open all hook and loop closures and slide the hand into the brace.
2. Make sure the brace is positioned below the palmar crease and the hand is seated comfortably in the palmar stay. Readjust and bend palmar stay as necessary.
3. Ensure the dorsal stay is centered over the top of the wrist. Readjust by removing the dorsal stay and aligning it over the center of the wrist.
4. Secure the hand closure between the thumb and index finger to the desired level of comfort.
5. Pull the lace tab away from the brace to fully tighten, then wrap over and around the wrist to secure.




NOTE: The brace should fit snug, but remain comfortable.

CLEANING INSTRUCTIONS

Hand wash with mild detergent in lukewarm water. Air dry.

WARRANTY

DeRoyal products are warranted for ninety (90) s from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**



WARRIOR® WRIST LACER

	NOT MADE WITH NATURAL RUBBER LATEX
	SINGLE PATIENT USE ONLY
	NON-STERILE
RX ONLY	Federal U.S.A. law restricts this device to sale or use by or on the order of a physician (or properly licensed practitioner).

IMPORTANT INFORMATION

Please read all instructions, warnings and precautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE

The DeRoyal® Warrior® Wrist Lacer brace is intended to provide compression, support and alignment to the wrist

INDICATIONS




Carpal tunnel syndrome; Tendonitis; Sprains and strains of the wrist; Overuse injuries

CAUTION

- This product is to be fitted initially by a physician or properly licensed practitioner who is familiar with the purpose for which it is intended. The physician or practitioner is responsible for providing wearing instructions and precautions to other healthcare practitioners or healthcare providers involved in the patient's care as well as the patient. Consult your physician or properly licensed practitioner immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this product. Discontinue use if pressure ulcers develop. Follow the instructions of your physician or properly licensed practitioner for duration of use. These instructions do not supersede hospital protocol and/or direct orders of patient's physician.
- Discontinue use of the device and consult physician or properly licensed practitioner if it becomes damaged in any way, limiting its ability to perform as intended.
- Inspect brace for damaged or missing components prior to patient application.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.**

INSTRUCTIONS FOR USE

1. Open all hook and loop closures and slide the hand into the brace.
2. Make sure the brace is positioned below the palmar crease and the hand is seated comfortably in the palmar stay. Readjust and bend palmar stay as necessary.
3. Ensure the dorsal stay is centered over the top of the wrist. Readjust by removing the dorsal stay and aligning it over the center of the wrist.
4. Secure the hand closure between the thumb and index finger to the desired level of comfort.
5. Pull the lace tab away from the brace to fully tighten, then wrap over and around the wrist to secure.



NOTE: The brace should fit snug, but remain comfortable.

CLEANING INSTRUCTIONS

Hand wash with mild detergent in lukewarm water. Air dry.

WARRANTY

DeRoyal products are warranted for ninety (90) s from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**



Manufacturer

DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 U.S.A.
888.938.7828
865.938.7828

Part#0-1940 Revised 1/17

©2017 DeRoyal Industries, Inc.
All Rights Reserved. Warrior, DeRoyal, the DeRoyal
logo, and Improving Care. Improving Business. are
trademarks or registered trademarks of
DeRoyal Industries, Inc.

Manufacturer

DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 U.S.A.
888.938.7828
865.938.7828

Part#0-1940 Revised 1/17

©2017 DeRoyal Industries, Inc.
All Rights Reserved. Warrior, DeRoyal, the DeRoyal
logo, and Improving Care. Improving Business. are
trademarks or registered trademarks of
DeRoyal Industries, Inc.