

Crystalline™ ECO

SUCTION SYSTEM



Environmentally Conscious Rigid Canister System

- Patented Design
- Hydrophobic Filter Technology
- Identification Icons on Lid
- Disposable System

Reduces Risk of Failure and Maintains Integrity for an **Average of 11x Longer*** when exposed to lipids in bodily fluids.




Competitor
Canister A



Competitor
Canister B




DeRoyal's ECO
Canister




Traditional Rigid Design

Lipids were able to penetrate the pores and expand the canister walls.

When interior surface of walls were under tension, this led to increased risk of failure.





Lid Enhancements ECO Rigid Canister System

Ordering Information

Product #	Description	Qty
71-0800	800cc Rigid Canister	50/Cs
71-8002	1200cc Rigid Canister, 2 elbows	40/Cs
71-8005	1200cc Rigid Canister w/DISS Mount, 1 elbow	40/Cs
71-8003	1500cc Rigid Canister, 2 elbows	40/Cs
71-8004	2500cc Rigid Canister, 2 elbows	32/Cs





DeRoyal Industries, Inc.
200 DeBusk Lane, Powell, TN 37849 USA
888.938.7828 | deroyal.com
Reprint# LIT0-2588 | Rev. 1/2023

The DeRoyal logo and Crystalline are registered trademarks of DeRoyal Industries, Inc.