



# Suction Tubing, Handles and Connectors

---

SOLUTIONS

# Suction Tubing and Connectors

## Suction Tubing with Connector

- Molded female connectors provides secure connection
- Male connector included for coupling additional tubes
- Minimal coil memory allows for easy handling and drapability
- Color coded connectors for size identification
- Available in economical non-sterile 100ft. rolls
- 100 ft rolls with cut-to-size female connectors every 6 ft
- 16 rigid male connectors included
- Not made with natural rubber latex

### $\frac{3}{16}$ " Tubing with Green Female Connectors



Product #	Description	Quantity
71-C501*	$\frac{3}{16}$ " x 1.5' (5 mm x .46 m)	100/Cs
71-C506*	$\frac{3}{16}$ " x 6' (5 mm x 1.8 m)	50/Cs
71-C510*	$\frac{3}{16}$ " x 10' (5 mm x 3 m)	30/Cs
71-C512*	$\frac{3}{16}$ " x 12' (5 mm x 3.7 m)	30/Cs
71-C520*	$\frac{3}{16}$ " x 20' (5 mm x 6.1 m)	25/Cs
71-C5100NS*	$\frac{3}{16}$ " x 100' (5 mm x 30 m)	1 Rl/Bx

### $\frac{1}{4}$ " Tubing with Light Blue Female Connectors





Product #	Description	Quantity
71-C606*	$\frac{1}{4}$ " x 6' (6 mm x 1.8 m)	45/Cs
71-C610*	$\frac{1}{4}$ " x 10' (6 mm x 3.0 m)	25/Cs
71-C612*	$\frac{1}{4}$ " x 12' (6 mm x 3.7 m)	25/Cs
71-C620*	$\frac{1}{4}$ " x 20' (6 mm x 6.1 m)	20/Cs
71-C6100NS*	$\frac{1}{4}$ " x 100' (6 mm x 30 m)	1 Rl/Bx

### $\frac{9}{32}$ " Tubing with Dark Blue Female Connectors




Product #	Description	Quantity
71-C706*	$\frac{9}{32}$ " x 6' (7 mm x 1.8 m)	40/Cs
71-C710*	$\frac{9}{32}$ " x 10' (7 mm x 3.0 m)	20/Cs
71-C712*	$\frac{9}{32}$ " x 12' (7 mm x 3.7 m)	20/Cs
71-C720*	$\frac{9}{32}$ " x 20' (7 mm x 6.1 m)	15/Cs
71-C7100NS*	$\frac{9}{32}$ " x 100' (7 mm x 30 m)	1 Rl/Bx




 **WARNING**- This product can expose you to chemicals including di(2-ethylhexyl)phthalate, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

# Suction Tubing and Connectors


## Yankauer Suction Handles



71-D880



71-D886



71-D887



71-D888

- Constructed of a shatter-resistant plastic
- Smooth inner surface for maximum flow
- Slip-resistant and lightweight, provides optimal control and maneuverability
- Translucent material allows for visual monitoring and inspection during use
- Available in three tip styles: bulbous, tapered bulbous and straight
- Bulbous tip is rounded to help minimize tissue grab and/or trauma
- Tapered bulbous tip helps minimize tissue grab and/or trauma, while allowing for better visibility
- Straight tip provides less bulk, allowing for better visibility
- Each tip has side eyes to help reduce tissue grab
- Available with control vent for intermittent or continuous suction control
- Available with pre-attached tubing
- Not made from natural rubber latex

Product #	Description	Quantity
71-D886	Bulbous tip, non-vented	50/Cs
71-D886V	Bulbous tip, vented	50/Cs
71-D880	Tapered bulbous tip, non-vented	50/Cs
71-D882	Tapered bulbous tip, vented	50/Cs
71-D887	Straight tip, non-vented	50/Cs
71-D887V	Straight tip, vented	50/Cs
71-D888	Bulbous Tip with On/Off Switch	50/Cs
71-D883A*	Tapered bulbous tip, non-vented with ¼" (6mm) x 6' (1.8m) tubing	50/Cs
71-D884*	Tapered bulbous tip, non-vented with ¼" (6mm) x 9' (2.7m) tubing	25/Cs
71-D885*	Tapered bulbous tip, non-vented with ¼" (6mm) x 12' (3.7m) tubing	25/Cs



71-D883A

# Suction Tubing and Connectors

## Flexible Suction Handles



71-D860

- Constructed of a malleable non-breakable vinyl with smooth inner lumen
- Highly flexible for precise suctioning
- Transparent for visual monitoring and inspection during use
- Available in three tip styles for different volumes of suction (fine, regular, high)
- Side eyes help to reduce the chance of tissue grab and/or trauma
- Vented option available with regular capacity
- Not made from natural rubber latex

Product #	Description	Quantity
71-D860	Fine capacity	50/Cs
71-D861	Regular capacity	50/Cs
71-D861V	Regular capacity, vented	50/Cs
71-D862	High capacity	50/Cs

## Poole Suction Handle




71-D8770

- Constructed from a shatter-resistant plastic
- Designed to remove pooled blood and fluids during surgical procedures
- Features three multipurpose positions that convert from sump action to pinpoint capabilities
- Outer sheath has 132 large eyes for rapid suctioning while reducing tissue grab
- Inner cannula provides pinpoint suctioning
- Universal connector for easy connection to a variety of tubing sizes
- Not made from natural rubber latex

Product #	Description	Quantity
71-D8770	Poole Suction Handle	50/Cs

## Frazier Suction Handles



71-D871

- Ideal for surgical procedures that require fine and delicate suctioning
- Malleable shaft can be bent to desired position without occlusion
- Obturator can be used to facilitate shaping of shaft and/or removal of obstructions that may become trapped during suctioning
- Obturator must be removed before suctioning begins
- Control vent allows for intermittent or continuous suction control
- Universal connector for easy connection to a variety of tubing sizes
- Not made from natural rubber latex

Product #	Description	Quantity
71-D871	8 French	50/Cs
71-D872	10 French	50/Cs
71-D873	12 French	50/Cs

# Suction Tubing and Connectors

## Connectors

### 5-in-1 Connector


- Compatible with a variety of medical tubing
- Fits tubing sizes  $\frac{3}{16}$  to  $\frac{7}{16}$  in. (5 to 11 mm)
- Available in sterile and non-sterile



Product #	Description	Quantity
71-515S	5-in-1 Straight Connector, sterile	1/Pch, 25 Pch/Bx, 8 Bx/Cs, 200/Ea
71-515N	5-in-1 Straight Connector, non-sterile	50/Bg, 1 Bg/Bx, 10 Bx/Cs, 500/Ea

### 6-in-1 Y Connector

- Compatible with a variety of medical tubing
- Fits tubing sizes  $\frac{1}{4}$  to  $\frac{5}{8}$  in. (6 to 15 mm)
- Available in sterile and non-sterile



Product #	Description	Quantity
71-616S	6-in-1 Y Connector, sterile	1/Pch, 15 Pch/Bx, 8 Bx/Cs, 120/Ea
71-616N	6-in-1 Y Connector, non-sterile	15/Bg, 1 Bg/Bx, 10 Bx/Cs, 150/Ea

### Sims Connector

- Compatible with a variety of medical tubing
- Fits tubing sizes  $\frac{1}{8}$  to  $\frac{1}{4}$  in. (3 to 6 mm)
- Available in sterile and non-sterile



Product #	Description	Quantity
71-618S	Sims Connector, sterile	1/Pch, 50 Pch/Bx, 8 Bx/Cs, 400/Ea
71-618N	Sims Connector, non-sterile	50/Bg, 1 Bg/Bx, 10 Bx/Cs, 500/Ea

### T Connector

- Compatible with a variety of medical tubing
- Fits tubing sizes  $\frac{3}{16}$  to  $\frac{7}{16}$  in. (5 to 11 mm)
- Available in sterile and non-sterile



Product #	Description	Quantity
71-619S	T Connector, sterile	1/Pch, 15 Pch/Bx, 8 Bx/Cs, 120/Ea
71-619N	T Connector, non-sterile	15/Bg, 1 Bg/Bx, 10 Bx/Cs, 150/Ea

### Elbow Right Angle Connector

- Connect suction tubing to ports at a 90° angle
- Available in sterile and non-sterile



Product #	Description	Quantity
71-1118	Elbow Right Angle Connector	100/Cs



DeRoyal Industries, Inc.  
200 DeBusk Lane, Powell, TN 37849 USA  
888.938.7828 | DeRoyal.com