

Prospera Spectrum NPWT System

The Prospera Spectrum NPWT System redefines the gold standard in NPWT performance. The system integrates Prospera Smart Propel™ Technology with an intuitive interface for use in acute and post acute settings. It is indicated for use on a variety of wound types, including acute and chronic wounds.


- » Multi-lumen Prospera Smart Propel Technology
- » GPS / RFID Location Tracking Technology
- » Pressure ranging from 20 to 200 mmHg¹
- » Clear and intuitive user interface
- » Continuous & Variable Pressure Modes
- » Audible and visual alarms
- » Long battery life (8+ hours)¹
- » Quiet operation




Redefining the GOLD Standard in NPWT Performance

- » Prospera Smart Propel Technology enables efficient fluid removal, supporting best in class wound pressure management¹
- » Prospera Spectrum NPWT system’s intuitive user interface supports ease of use in both acute and post acute settings
- » Compatible with DeRoyal’s Continuum® Inventory Management System, facilitating efficient care from hospital to home

Average Canister Fill Efficiency over 12 Days¹



Wound Pressure Variability 125 mmHg² Normal Viscosity



Scan QR Code for the Prospera Spectrum NPWT System Overview

- » The height of each plot shows the range of pressures experienced at the wound bed. The width shows the amount of measurements recorded at that pressure value
- » The Prospera Spectrum NPWT System most consistently delivers therapy within the range of the prescribed pressure¹
- » This trend can be observed across varying conditions such as: low or high exudate rate, thin or viscous exudate, and presence or absence of a leak²



Product No.	Description	Packaging
NP-7000	Prospera Spectrum NPWT System Unit, Power Supply, Accessories	1 Ea
NP-7000-GPS	Prospera Spectrum NPWT System w/ GPS Unit, Power Supply, Accessories	1 Ea
NP-7300	Prospera Spectrum 300cc Canister w/ solidifier	10 Ea/Cs
NP-7500	Prospera Spectrum 500cc Canister w/ solidifier	10 Ea/Cs
NP-7202	Prospera Spectrum Charging Cord	1 Ea
NP-7205	Prospera® NPWT Unit IV / Bed Rail Clamp	1 Ea
NP-6006U	Small Kit: Black Foam (10cm x 8cm x 3cm), Prospera Smart Propel™ Dome, 1 Transparent Acrylic Drape (12" x 12"), Ruler	10 Ea/Cs
NP-6007U	Medium Kit: Black Foam (18cm x 13cm x 3cm), Prospera Smart Propel Dome, 2 Transparent Acrylic Drape (12" x 12"), Ruler	10 Ea/Cs
NP-6008U	Large Kit: Black Foam (26cm x 15cm x 3cm), Prospera Smart Propel Dome, 3 Transparent Acrylic Drape (12" x 12"), Ruler	10 Ea/Cs
NP-0080U	X-Large Kit: Black Foam (45cm x 30cm x 3cm), Prospera Smart Propel Dome, 4 Transparent Acrylic Drape (12" x 12"), Ruler	5 Ea/Cs

To learn more about DeRoyal's products or to place an order, Call 1-888-938-7828 or visit DeRoyal.com.

¹ Data on File

² Data on File

DeRoyal, the DeRoyal logo, Prospera, Prospera Spectrum, and Prospera Smart Propel are trademarks or registered trademarks of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | 1.888.938.7828 | deroyal.com

© 2026, DeRoyal. All rights reserved. | Reprint # 0-2643 | Rev. 3/2026