

Introducing DeRoyal's Prone Patient Positioning Kit

Provides a safe universal solution for prone positioning.

Head-to-toe care, support, stability, and comfort for patient protection when placed in a prone position. This kit is intended to be used on patients who may potentially benefit from the prone position, such as patients with acute respiratory distress syndrome, acute hypoxemia, and non-cardiac pulmonary edema.

- Helps reduce risk of hospital acquired pressure injuries
- Aids in safe consistent alignment of patient
- Can be used for ICU prone patient positioning
- Not made with natural rubber latex



 DeRoyal®

Prone Patient Positioning Kit



Kit contains:



Head Positioner

10" x 10.5" x 2.5" Soft Head Positioner with "T" cutout to allow for secure endotracheal tube placement



Standard Foam Pads

8" x 11" x 1" Standard high-density foam with soft top layer for increase patient comfort




Large Foam Pad

15" x 11" x 1" Large high-density foam with soft top layer for increase patient comfort



Shin Positioners

9" x 6" Semicircular high-density shin positioner



Ordering Information

Product #	Quantity
M60-360	5 Kits/Cs



Improving Care. Improving Business.®